

DISCOVER THE POWER OF CONNECTION

MINIMED 640G®





INSULIN PUMP THERAPY CHANGING LIVES TODAY



Has your child just been diagnosed with insulin dependent diabetes? Or perhaps they've been on multiple daily injection therapy for years?

Whatever the situation, Medtronic is here to support your family in better managing diabetes with insulin pump therapy.

An insulin pump is a discreet electronic device, about the size of a small mobile phone. It closely mimics the way a healthy pancreas functions, replacing the need for frequent injections by delivering precise doses of rapid-acting insulin 24 hours a day to closely match their body's daily needs.

HOW DOES THE PUMP WORK?

The pump provides insulin to your body at two rates:



BASAL RATE

Small amounts of insulin released continuously throughout the day to mimic the background insulin production of the pancreas.



BOLUS RATE

Additional insulin can be delivered "on demand" to match the food you are going to eat or correct high blood sugar.

COMPONENTS OF INSULIN PUMP THERAPY



1. INSULIN PUMP

A small, durable electronic device to program your insulin & show you how you are tracking. The device includes a reservoir compartment.





2. INFUSION SET

An infusion set includes a thin cannula that goes from the reservoir to the infusion site on your body. The cannula is inserted into the site - similar to where you would give insulin injections. The infusion set is changed every two to three days.



3. RESERVOIR

A plastic cartridge that holds he insulin and is locked into the nsulin pump. A reservoir can hold up to 300 units of insulin and is changed every two to three days.

INSULIN PUMP THERAPY THE BENEFITS



With an insulin pump like the MiniMed[®] 640G, you can better manage the need for insulin dose adjustment (particularly for meals, exercise, on sick days and overnight) helping your child achieve better glucose control.

It is easier to do things like eating out, sleeping-in on the weekends, playing sports, travelling and even swimming as the MiniMed[®] 640G is waterproof up to 3.6m for 24 hours.

FEATURE	PUMP	MDI
Fewer Daily Injections	1	×
24/7 Insulin Infusion	1	×
High and Low Alerts & Alarms	1	×
Insulin Suspension to Avoid Hypos	1	×
Bolus Calculator	1	×
Ongoing Insulin Tracking to Avoid Stacking	1	×
Visualise your Glucose Data	1	×

FEWER HYPOGLYCAEMIC EVENTS¹



EASIER DOSING. ACCURATE DOSING. NO NEEDLES

GREATER FLEXIBILITY IN WHEN TO EAT AND HOW TO EXERCISE

THE MINIMED[®] 640G HOW DOES IT WORK?



PROTECTING YOU TO STAY IN TARGET

The MiniMed[®] 640G System³ keeps you **in target range** with custom basal profiles and an in-built bolus wizard.

PROTECTING YOU AGAINST HYPOS

The MiniMed[®] 640G System³ is the only system that offers **advanced protection from lows** with predictive low glucose management technology, especially throughout the night.

PROTECTING YOUR QUALITY OF LIFE

The MiniMed[®] 640G System³ **discreetly fits** in with your lifestyle with its waterproof design and remote bolus for snacks.

HOW SMARTGUARD[™] WORKS[^]



^SmartGuard Technology is activated when the MiniMed 640G is paired with continuous glucose monitoring *The dynamic suspend feature is based on certain criteria: sensor glucose must be within 3.9mmol/L of the low limit and predicted to be 1.1mmol/L above the low limit within 30 minutes AND the pump must not be in the refractory period. ** The dynamic resume feature is based on certain preset criteria: sensor glucose must be 1.1mmol/L above the preset low limit and predicted to be 2.2mmol/L above within 30 minutes AND insulin must have been suspended for at least 30 minutes.

INSULIN DELIVERY FOR EVERY BODY



The infusion set consists of tubing that carries insulin from the pump to you.

On one end of the tubing is the reservoir connector that is attached to the reservoir which holds the insulin. On the other end is the insertion site section that attaches to you.

Make sure to change your infusion sets every 2-3 days to ensure continuous insulin delivery.

Medtronic offers four (4) different types of infusion sets for you to choose from:

Sure-T[®] infusion set

A simple steel needle infusion set with 90° insertion which is suitable for those who have had reactions to plastic cannula or a history of bent cannula. Also suitable for pregnant women (up to the 2nd trimester) and small children.





Quick-set[™] infusion set

Our most popular plastic cannula infusion set with 90° insertion. Suitable for those who are less active, have limited dexterity and for alternative site usage.

Silhouette[™] infusion set

A flexible insertion angle plastic cannula infusion set with 30-45° insertion. Suitable for those who are active, lean or muscular. Also suitable for pregnant women (up to the 2nd trimester).





mio[®] infusion set A convenient 'all in one' plastic cannula infusion set with 90° insertion. Suitable for those who are less active, and for alternative site usage.

GET THE FULL PICTURE WITH CARELINK



Our three-step approach to diabetes management is fully integrated. Your Medtronic pump will wirelessly connect to your CONTOUR® NEXT LINK 2.4 blood glucose meter, meaning all data is stored conveniently on your pump.



Our Carelink[®] Personal Software helps you visualise your diabetes information with easy to read charts, statistics and events that can help you identify and understand patterns and trends. With one simple upload, you and your healthcare team can make adjustments to your therapy or remotely.

Visit: carelink.medtronic.eu for more information



Example Carelink Report

REAL STORIES FROM REAL PEOPLE[‡]



At just 29 months old, Shiloh was diagnosed with Type 1 diabetes. An insulin pump has changed the way her parents manage her condition

"To get a pump was the greatest thing we've ever done really... she can have the freedom of a normal child!"

"The benefits I see with the pump is that she can do anything that a normal child can do - dancing, swimming, whatever she wants to."

"My message to other parents is don't be afraid of it - you can live life first and have diabetes follow." ~ Shiloh's parents

Check out more patient stories by visiting: medtronic-diabetes.com.au



By choosing a MiniMed[™] 640G Insulin Pump you have joined the global Medtronic Diabetes Family.

For a full range of accessories to style your pump visit : eshop.medtronic-diabetes.com.au

REFERENCES

1. Bergenstal RM1, Tamborlane WV, Ahmann A, Buse JB, Dailey G, Davis SN, Joyce C, Perkins BA, Welsh JB, Willi SM, Wood MA; STAR 3 Study Group. Sensor-augmented pump therapy for A1C reduction (STAR 3) study: results from the 6-month continuation phase. Diabetes Care. 2011 Nov;34(11):2403-5. doi: 10.2337/dc11-1248. Epub 2011 Sep 20.

2. Assumes four injections per day for 30 days and one infusion set change every three days.

3. The system is inclusive of the MiniMed 640G Insulin Pump, Continuous Glucose Monitoring System (Guardian 2 Link transmitter ans Enlite sensor), Contour Next LINK 2.4 blood glucose meter and Medtronic CareLink therapy management software.

MiniMed, Guardian & CareLink are trademarks of Medtronic MiniMed, Inc. © 2017 Medtronic MiniMed, Inc. All Rights Reserved.

‡ This patient testimonial relates to an account of an individual's response to treatment. The account is genuine, typical and documented. The response other persons have to treatment could be different. Please consult with you healthcare professional for information on whether the treatment is appropriate for you.

Medtronic

Medtronic Diabetes Australia 2 Alma Road Macquarie Park, New South Wales, 2113 Australia

24 Hour Global Helpline 1800 777 808

medtronic-diabetes.com.au