

# Which infusion set is right for you?

## STEEL

- Reactions to plastic cannula
  - History of bent cannula
- Pregnant (up to 2nd trimester)
  - Small child

90 degree insertion

**Sure-T® infusion set**  
Simple steel needle infusion set



**Low/Average BMI**

**Average/High BMI**

**6 mm**

**8 mm**

**10 mm**

Wear insulin pump close to infusion site  
**60 cm**

Wear insulin pump close to infusion site  
**60 cm**

**60 cm**

Wear insulin pump away from infusion site  
**80 cm**

Wear insulin pump away from infusion site  
**80 cm**

## PLASTIC CANNULA

- Less active
- Alternative site usage
- Limited dexterity

90 degree insertion

**Quick-set™ infusion set**  
Our most popular infusion set



**Low/Average BMI**

**Average/High BMI**

**6 mm**

**9 mm**

Wear insulin pump close to infusion site  
**60 cm**

Wear insulin pump close to infusion site  
**60 cm**

Wear insulin pump away from infusion site  
**80 or 110 cm**

Wear insulin pump away from infusion site  
**80 or 110 cm**

- Less active
- Alternative site usage

90 degree insertion

**mio® infusion set**  
Convenient "all in one" infusion set



**Low/Average BMI**

**Average/High BMI**

**6 mm**

**9 mm**

Wear insulin pump close to infusion site  
**45 or 60 cm**

**80 cm**

Wear insulin pump away from infusion site  
**80 cm**

- Active
- Lean
- Muscular
- Pregnant (from 2nd trimester)

30-45 degree insertion

**Silhouette™ infusion set**  
Flexible insertion angle infusion set



**Low/Average BMI**

**Average/High BMI**

**13 mm**

**17 mm**

Wear insulin pump close to infusion site  
**45 or 60 cm**

Wear insulin pump close to infusion site  
**60 cm**

Wear insulin pump away from infusion site  
**80 or 110 cm**

Wear insulin pump away from infusion site  
**80 or 110 cm**